

HIGH VOLUME SHIPPER / SMARTPAK II®

ADVANCED THERAPY SHIPPER™

SIDE 1

BATTERY RESET, REPLACEMENT & RECHARGE INSTRUCTIONS


STEP 1

Open the Smartpak II[®] Housing located on the lid of the high volume shipper by unlatching the black latch and pivoting the lid upwards.



STEP 2

Partially lift the Smartpak II[®] from the foam cradle. Disconnect the thermocouple from the mini-USB port.



STEP 3

Lift the Smartpak II[®] up out of the container.





TO REASSEMBLE UNIT FOLLOW THIS PROCEDURE IN REVERSE ORDER

IF THE DRY SHIPPER DOES NOT EMIT VAPOR WHEN PLUG IS REMOVED PLEASE CONTACT CUSTOMER SERVICE IMMEDIATELY, BY CALLING (949) 470-2305 OR EMAIL CS@CRYOPORT.COM




ADVANCED THERAPY SHIPPER™


BATTERY RESET, REPLACEMENT & RECHARGE INSTRUCTIONS


```
TO RESET
```

Default configuration of the Smartpak II® is to record and transmit sensor data at 5-minute and 60-minute intervals, respectively when it is connected to a stable cellular network. The Smartpak II® device can be effectively reset per the instructions below.

Continuously tap and release the touch sensitive button at the frequency of one tap per 3 seconds until the Network LED light turns on — regardless of what the color is (Red or Green). It should take approximately 5 taps for the Network LED light to turn on.







SIDE 2

Wait to see a solid Green LED light. After 5 minutes, refresh your Live View[™] link. The queued-up data should fully be displayed up to the timestamp you executed the step above.

PLEASE NOTE: The Network LED light may be "blinking" RED. The blinking could last up to one minute (depending on the cellular coverage condition) before it changes to a solid Green.

