

ADVANCED THERAPY SHIPPERS™

EXPANDED TEMPERATURE SOLUTIONS



CRYOPORT. CERTIFIED. COOL.™

Cryoport is proud to introduce our new C3™ Advanced Therapy Shippers™. Our Advanced Therapy Shippers™ are dedicated exclusively to cell and gene supply chain management for human clinical and commercial therapies. For valuable commodities that need to be kept at Refrigerated or Controlled-Room Temperatures, the C3™ Solution is the most advanced temperature controlled shipper in the life science market. The validated C3™ shipper is comprised of the Smartpak II® Condition Monitoring System, Cryoport® Logistics Management Platform and 24/7/365 monitoring and customer support. The industry-leading Smartpak II® and powerful Cryoport® work in conjunction to provide visibility of the location and the key aspects of your critical shipment. Cryoport's unparalleled technology and round-the-clock customer support allows for risk mitigation, proactive monitoring and intervention capabilities.

Specs	Refrigerated	CRT
Temperature Range	2 – 8°C	15 – 25°C
Duration – Summer Profile	96 hours	96 hours
Duration – Winter Profile	96 hours	72 hours
Volume	12 L / 729 ci	
Tare Weight	22 lbs / 9.98 kg	
Payload Size Imperial	8.9 x 8.9 x 8.9 inches	
Payload Size Metric	226 x 226 x 226 mm	
Outer Dimension Imperial	16 x 14 x 15 inches	
Outer Dimension Metric	407 x 356 x 381 mm	



The Smartpak II® Condition Monitoring System is seamlessly integrated inside the C3™ shipper and monitors: location, temperature, shock, orientation, pressure, humidity, anti-tamper and battery life.

KEY FEATURES:

Packaging

- New and recyclable outer packaging for every shipment to provide a clean and consistent appearance
- Simple pack-out with unique flap design for easy labeling identification
- Stationary temperature probe to maintain consistent location for measurement
- Serialized phase change panels and cooler are each tracked per shipment
- Phase change panels surround the payload providing greater temperature performance and overall payload protection
- C3™ shippers consist of durable components that meet ISTA 3A transit criteria

Informatics

- The Cryoport Solution, comprised of the Smartpak II® and Cryoport® provide unsurpassed analytics for planning, enroute mitigation and reporting
- Cryoport® streamlines ordering, tracking, paperwork and communications through a single web-based portal
- Chain of Compliance® processes provide full traceability and transparency of the commodity history, historical use, performance, handling, cleaning and requalification of the equipment
- Complete temperature record history available in 21 CFR Part 11 compliant database

Logistics

- Cryoport's unparalleled expertise in cold chain logistics solutions and our exclusive focus on medical transport
- 24/7/365 customer service support to track your shipment and mitigate risk with automated escalation
- Flexibility on freight carrier options
- Carrier optimization and route planning aimed at balancing risk and performance cost
- Cryoport provides complete chain of custody, chain of condition and Chain of Compliance® with every shipment



SCIENCE. LOGISTICS. CERTAINTY.

U.S. • The Netherlands • Singapore
info@cryoport.com | +1 949.232.1900 | cryoport.com